

A Flow Cytometric, Whole Blood Dendritic Cell Immune Function Assay

Cross-reference to Related Applications

This application is a continuation-in-part of
5 commonly owned and copending U.S. Patent application
serial no. 09/158,406, filed September 22, 1998, the
disclosure of which is incorporated herein by reference
in its entirety.

Field of the Invention

10 This invention relates to assays of blood
cell function, and in particular to assays of dendritic
cell function in whole blood.

Background of the Invention

Dendritic cells (DCs), first identified a
15 quarter century ago by a characteristic "dendritic"
morphology observable in peripheral lymphoid tissues,
Steinman *et al.*, *J. Exp. Med.* 137:1142-1162 (1973), are
now known to be a morphologically-diverse and widely-
distributed cell population. Today, these diverse
20 cells are collectively distinguished by a common

function: dendritic cells are the most potent antigen-presenting cells (APCs) of the mammalian immune system, and alone among the various antigen-presenting cells appear capable of triggering a primary T lymphocyte
5 response.

This singular ability to prime a T cell-mediated immune response — combined with a potent ability to present antigen to activated T cells — has implicated dendritic cells as potential reagents for
10 immune-based therapies, as well as likely targets for therapeutic intervention in the treatment of various immune-mediated disorders.

For example, WO 97/24438 describes compositions and methods for co-culturing dendritic
15 cells with T lymphocytes and protein antigen *in vitro*, thus driving the *ex vivo* antigen-specific activation of T cells. The activated T cells are then administered autologously to effect an antigen-specific immune response *in vivo*. Similarly, WO 97/29183 describes a
20 method of activating T cells *in vitro* by contacting the T lymphocytes with DC that directly express an antigenic protein from a recombinant construct. Again, the activated T cells are intended for autologous infusion. Specific application of DC-driven *ex vivo* T
25 cell activation to the treatment of prostate cancer is described and claimed in U.S. Patent No. 5,788,963. In yet another approach, Nemazee, U.S. Patent No. 5,698,679, describes and claims immunoglobulin fusion proteins that deliver antigenic peptides to targeted
30 antigen presenting cells (APCs), including dendritic cells, *in vivo*.

Dendritic cells have also been implicated as important in the pathogenesis and pathophysiology of

AIDS. One type of DC, the Langerhans cells (LC), is generally believed to be the initial cell type infected with HIV following mucosal exposure to virus. DC are believed to act not only during the initial phase of
5 HIV disease, but also during the chronic phase, facilitating infection and depletion of T lymphocytes. Zoetewij et al., *J Biomed Sci* 5(4):253-259 (1998). DCs in lymphoid mucosa may represent a key reservoir of viral nucleic acid and virions throughout the course of
10 the disease. Grouard et al., *Curr. Opin. Immunol.* 9(4):563-567 (1997); Weissman et al., *Clin. Microbiol. Rev.* 1997 10(2):358-367 (1997). *In vitro* methods for screening pharmaceutical candidates for agents that abrogate HIV infection of DC are described and claimed
15 in Steinman et al., U.S. Patent No. 5,627,025.

Yet despite their importance to the normal mammalian immune response and in immunopathology, DCs have been difficult to study, and particularly difficult to study in their native milieu.

20 The difficulty stems in part from the rarity of dendritic cells. Although widely distributed, DC are sparse, even in lymphoid tissues, and represent no more than about 0.3% - 0.5% of nucleated cells in human peripheral blood.

25 A further difficulty arises from the absence of DC-specific cell surface markers that would readily permit the positive immunoselection of DCs from mixed populations of cells.

Extensive efforts to identify surface markers
30 that define DCs have been only partially successful. As a result, DCs are presently identified by multiple-marker panels, with identification based primarily on the absence of staining with markers for other lineages

(i.e., as lin^- cells). The result is that typical DC immunopurification protocols require at least one immunodepletion step, eliminating cells of various nondendritic blood lineages — lymphocyte, monocyte, granulocyte, and NK lineages, e.g. — coupled with at least one immunoenrichment step. The immunoenrichment step may, for example, include selection for CD4^+ cells (Blood Dendritic Cell Isolation Kit, Miltenyi Biotec #468-01, Auburn, CA), or, in the alternative or in addition, selection for HLA-DR expression, Ghanekar et al., *J. Immunol.* 157:4028-4036 (1996).

These serial manipulations, however, may substantially alter the DC cell phenotype from that present *in vivo*. For example, $\text{lin}^- \text{HLA-DR}^+ \text{CD123}^+$ dendritic cells in fresh preparations of tonsillar mononuclear cells express low levels of the T cell costimulatory molecules CD80 (B7.1), CD86 (B7.2), and HLA-DQ. Even an overnight culture of these cells in the absence of added cytokines is sufficient to induce the mature DC phenotype with upregulation of CD86, CD80, HLA-DQ and HLA-DR. Olweus et al., *Proc. Natl. Acad. Sci. USA* 94(23): 12551-12556 (1997). Longer term culture of CD34^+ dendritic cell precursors in the presence of cytokines effects substantial phenotypic changes. Caux et al., *J. Exp. Med.* 184:695, 1996; Olweus et al., *Proc. Natl. Acad. Sci. USA* 94(23):12551-12556 (1997).

Thus, there exists a need in the art for methods of assaying dendritic cells without prior immunopurification or *in vitro* culture.

The paucity of DC-specific cell surface markers further suggests that surface immunophenotypic markers may only incompletely distinguish dendritic

cell subsets that are, nonetheless, functionally distinct. For example, peripheral blood dendritic cells have been shown to fall into two subsets distinguishable by the divergent expression of CD11c and CD123: one subset is CD11c⁺CD123^{low}, the other CD11c⁻CD123⁺. Olweus et al., *Proc. Natl. Acad. Sci. USA* 94(23): 12551-12556 (1997). Yet the critical and disparate roles that dendritic cells play in the immune system would argue that these two subsets each likely encompasses a variety of cell types with disparate functional activity.

There thus exists a need in the art for methods of distinguishing dendritic cell subsets using phenotypic criteria other than, or in addition to, expression of cell-surface markers. There further exists a need for methods of subsetting DC based on criteria that may be related more directly to DC function.

Recently, several groups have reported that intracellular staining of cells using cytokine-specific antibodies permits the flow cytometric analysis of cytokine expression in highly purified blood cell lineages, including purified dendritic cells. Picker et al., *Blood* 86(4):1408-1419 (1995); Waldrop et al., *J. Clin. Invest.* 99:1739-1750 (1997); Ghanekar et al., *J. Immunol.* 157:4028-4036 (1996); de Saint-Vis et al., *J. Immunol.* 160:1666-1676 (1998). More recently, Suni et al., *J. Immunol.* 212:89-98 (1998) described an assay for concurrent expression of intracellular cytokines and cell surface proteins in antigen-stimulated T lymphocytes without prior T cell purification. Similar assays are described and claimed in co-owned and

copending U.S. patent application nos. 08/760,447 and 08/803,702.

There exists a need in the art for a method that would adapt intracellular cytokine assays to the measurement of cytokine production by unpurified DC cells in whole blood.

Summary of the Invention

The present invention solves these and other problems in the art by providing, in a first aspect, a flow cytometric method for measuring dendritic cell function in whole blood, comprising the steps of: (a) contacting a whole blood sample with a dendritic cell activator; (b) contacting the sample with a plurality of dendritic cell-distinguishing antibodies and at least one cytokine-specific antibody; and then (c) flow cytometrically assaying the sample for the binding of cytokine-specific antibody by at least one distinguishable DC subset.

In preferred embodiments, activation is performed in the presence of an inhibitor of protein secretion, and following permeabilization of the cells cytokines are detected intracellularly. Thus, in a particularly preferred embodiment, the dendritic cell activator contacting step is performed in the presence of Brefeldin A, and the antibody contacting step itself comprises the steps, in order, of: (b1) adding a plurality of dendritic cell-distinguishing antibodies to the sample; (b2) lysing erythrocytes in the sample; (b3) permeabilizing nucleated cells in the sample; and then (b4) adding at least one cytokine-specific antibody to the sample.

The dendritic cell-distinguishing antibodies may include a plurality of non-DC lineage-specific antibodies. In such cases, it is particularly preferred that each of the non-DC lineage-specific
5 antibodies be conjugated to the identical fluorophore. When a plurality of non-DC lineage-specific antibodies is used, the dendritic cell-distinguishing antibodies further include an antibody specific for HLA-DR.

In a preferred embodiment, subsets of
10 dendritic cells are distinguishably labeled. In this embodiment, the dendritic cell-distinguishing antibodies include at least one antibody that binds differentially to the surface of the different dendritic cell subsets. Particularly preferred in this
15 embodiment is the use of antibody specific for CD11c or CD123.

The function of DCs can be characterized by their cytokine expression patterns and by the dynamic regulation of differentiation/activation markers (CMRF-
20 44, CMRF-56, CD83, CD25), of co-stimulatory molecules (CD40, CD80, CD86) and of class II major histocompatibility complexes (MHC class II). Thus, in a second aspect, the invention provides a flow cytometric method for measuring dendritic cell function
25 in whole blood, comprising: (a) contacting a whole blood sample with a dendritic cell activator; (b) adding to the sample a plurality of dendritic cell-distinguishing antibodies and at least one antibody specific for a dendritic cell surface marker indicative
30 of dendritic cell activation; and then (c) flow cytometrically assaying said sample for the binding of said antibody specific for the dendritic cell surface activation marker by at least one distinguishable DC subset. The dendritic cell surface activation marker

is usefully selected from the group consisting of differentiation/activation markers (CMRF-44, CMRF-56, CD83, CD25), of co-stimulatory molecules (CD40, CD80, CD86) and of class II major histocompatibility
5 complexes (MHC class II).

Brief Description of the Drawings

The above and other objects and advantages of the invention will be apparent upon consideration of the following detailed description, taken in
10 conjunction with the accompanying drawings, in which like reference characters refer to like parts throughout, and in which:

FIG. 1 is a flow chart schematizing the basic steps in a whole blood flow cytometric assay for
15 dendritic cell function, with LPS exemplified as the dendritic cell activator;

FIG. 2 presents a series of dot plots generated during the flow cytometric analysis of whole blood activated with LPS in the presence of
20 Brefeldin A. CD11c⁺ dendritic cells are painted green and CD11c⁻ dendritic cells are painted red; nondendritic cells appear gray. The colors are arbitrarily chosen, and bear no relationship to the fluorophores used for analysis;

25 FIG. 3 presents a series of dot plots generated during the flow cytometric analysis of whole blood activated with PMA+I in the presence of Brefeldin A. CD11c⁺ dendritic cells are painted green and CD11c⁻ dendritic cells are painted red; nondendritic
30 cells appear gray. The colors are arbitrarily chosen, and bear no relationship to the fluorophores used for analysis;

FIG. 4 presents a series of dot plots generated during the flow cytometric analysis of whole blood incubated in the presence of Brefeldin A in the absence of activator (resting control). CD11c⁺ dendritic cells are painted green and CD11c⁻ dendritic cells are painted red; nondendritic cells appear gray. The colors are arbitrarily chosen, and bear no relationship to the fluorophores used for analysis;

FIG. 5 presents the differential expression of TNF α , IL-8, CD80 and CD86 in CD11c⁺ dendritic cells from two donors, each activated alternatively with LPS or PMA+I;

FIG. 6 presents a series of histograms summarizing the effects of three different dendritic activators on the surface expression of the identified markers on peripheral blood dendritic cells in whole blood;

FIG. 7 shows a comparison of cytokine expression between monocytes (gray bars) and CD11c⁺ DCs (black bars) in activated whole blood, with FIG. 7A showing LPS + Brefeldin A-stimulated cells, and FIG. 7B showing PMA+I + Brefeldin A-stimulated cells; and

FIG. 8 shows kinetics of TNF α , IL-1 β , IL-6 and CD80 in LPS activated CD11c⁺ DCs. The time course of LPS incubation was 0 to 8 hours. Intracellular cytokine and CD80 expression is measured as PE mean fluorescence intensity (MFI).

Detailed Description of the Invention

In order that the invention herein described may be fully understood, the following detailed description is set forth. In the description, the following terms are employed:

By "whole blood" is intended a fluid blood sample as drawn from a mammal and substantially unfractionated thereafter. That is, if fractionation is performed subsequent to blood draw, the
5 fractionation has raised the percentage of dendritic cells to no more than about 5%, preferably no more than about 1 - 4%, most preferably no more than 1%, of total nucleated cells;

"Antibody" includes all products, derived or
10 derivable from antibodies or from antibody genes, that are useful as markers in the flow cytometric methods described herein. "Antibody" thus includes, *inter alia*, natural antibodies, antibody fragments, antibody derivatives, and genetically-engineered antibodies,
15 antibody fragments, and antibody derivatives;

"Dendritic cell-distinguishing antibody" includes any antibody that may be used, alone or in combination with other antibodies, to facilitate identification of dendritic cells, and thus includes
20 antibodies that are specific for epitopes displayed by non-DC lineages and further includes antibodies that bind to structures displayed by DC that prove useful for positive immuno-identification;

"Lineage negative", also abbreviated "lin-",
25 denominates the absence of cell surface markers known to be characteristic of non-dendritic lymphopoietic or hematopoietic cell lineages. By "absence" is intended a level of surface expression, as measured in an immunoassay, such as a flow cytometric assay, that is
30 not significantly different from background;

A "dendritic cell activator" is any substance that is capable of inducing or upregulating expression of cytokines, chemokines, or detectable cell surface proteins by dendritic cells;

All remaining terms have their usual meaning in the flow cytometric arts, as set forth, *inter alia*, in Ormerod (ed.), Flow Cytometry: A Practical Approach, Oxford Univ. Press (1997); Jaroszeski et al. (eds.),
5 Flow Cytometry Protocols, Methods in Molecular Biology No. 91, Humana Press (1997); and Practical Flow Cytometry, 3rd ed., Wiley-Liss (1995).

Dendritic cells (DC) capture, process and present antigen to naive and memory T cells, and thus
10 play a pivotal role in the mammalian immune response. An understanding of DC function is critical to any detailed understanding of mammalian immune function. Yet functional studies of dendritic cells have in the past been hampered by the functional diversity of the
15 cells that are collectively so denominated.

For example, studies of peripheral blood dendritic cells were for two decades conducted without awareness of the fact that peripheral blood dendritic cells fall into two mutually-exclusive subsets
20 distinguishable by cell surface immunophenotype. Thomas et al., *J. Immunol.* 153:4016 (1994); O'Doherty et al., *Immunology* 82:487-493 (1994); Olweus et al., *Proc. Natl. Acad. Sci. USA* 94(23): 12551-12556 (1997). Both subsets express high levels of HLA-DR and lack
25 markers characteristic of other lineages (CD3, CD14, CD19, CD20, CD16, CD56). The subsets are distinguished from one another by their divergent expression of CD11c and CD123: one subset is CD11c⁺CD123^{low}, the other CD11c⁻CD123⁺. O'Doherty et al., *Immunology* 82:487-493 (1994);
30 Olweus et al., *Proc. Natl. Acad. Sci. USA* 94(23): 12551-12556 (1997); Willmann et al., "Peripheral Blood

Dendritic Cells Revealed by Flow Cytometry" (Becton-Dickinson Application Note 3) (1998).

The two peripheral blood DC subsets that were first identified by fortuitous cell surface distinctions have now been shown to be functionally distinct. It is known, for example, that the CD11c⁺CD123^{low} DC subset proves more potent than the CD11c⁻CD123⁺ subset in stimulating T cells in a mixed lymphocyte reaction (MLR). And as shown newly herein, the CD11c⁺ subset alone responds to DC activators with upregulation of cytokine production and increased surface expression of T cell costimulatory molecules.

The two decades that intervene between the first identification of DC and the first demonstration that peripheral blood contains immunophenotypically and functionally discrete DC subsets speak to the insufficiency of surface phenotyping fully to capture the functional diversity of dendritic cells.

The present invention permits peripheral blood dendritic cells to be described and distinguished based upon differences in their functional responses to DC activators. The invention further permits these functional responses to be measured with minimal experimental intervention, precluding the known phenotypic plasticity of dendritic cells from confounding the results.

FIG. 1 schematizes the basic method of the present invention. A sample of whole blood is first incubated with a DC activator. LPS is exemplified in the figure.

Incubation with the dendritic cell activator serves to drive the differential phenotypic response of the various DC subsets present in the sample; measurement of these differences permits the

discrimination of DC subsets that might otherwise prove indistinguishable. Different activators produce different sets of responses, permitting still finer distinctions to be drawn. Although both FIG. 1 and the 5 experiments reported herein exemplify the invention using DC activators with broad and pleiotropic effects, such as LPS, activators with finer specificity will also prove useful.

Incubation with a dendritic cell activator is 10 particularly shown differentially to upregulate the production of various cytokines by peripheral blood dendritic cells which, absent stimulation, produce no detectable cytokines. By performing the activation 15 step in the presence of Brefeldin A ("BFA"), which disassembles the Golgi complex (Openshaw et al., J. Exp. Med. 182:1357 (1995); Chardin and McCormick, Cell 97:153 (1999)), inhibiting protein transport through the cellular secretion pathway, cytokine proteins 20 accumulate in a later step of the assay. Similar results would be obtained using equivalent inhibitors of secretion, such as monensin.

After incubation in the presence of activator and BFA, the surface of the cells is stained with 25 fluorophore-conjugated antibodies.

This surface staining step includes, as a first class of antibodies, a plurality of dendritic cell-distinguishing antibodies. A dendritic cell-distinguishing antibody is any antibody that may be 30 used, alone or in combination with other antibodies, to facilitate identification of dendritic cells. Thus, the antibodies used in this step may include (1) antibodies that preferentially bind non-dendritic

cells, and (2) antibodies that bind to dendritic cell surface structures useful in identifying DC.

As to the first such category, a cocktail of lineage-specific antibodies labeled with the identical fluorophore may advantageously be used. One such cocktail available commercially is the lin 1 FITC lineage cocktail from Becton Dickinson Immunocytometry Systems (BDIS, San Jose, CA, catalogue number 340546), which contains a mixture of antibodies specific for CD3, CD14, CD16, CD19, CD20, and CD56, each conjugated to fluorescein isothiocyanate (FITC). In combination, the antibodies in the cocktail stain lymphocytes, monocytes, eosinophils, and neutrophils, but not dendritic cells. The DC in the labeled sample thus assort into the FITC⁻ or FITC^{low} class. The lin 1 cocktail is particularly advantageous in that the concentration of antibodies and degree of conjugation have been titrated to provide equivalent intensity fluorescence signals from the cells of the various non-DC lineages that are bound by the antibodies.

As to the second category of dendritic cell-distinguishing antibodies there is, as yet, no cell surface marker that alone positively identifies dendritic cells. When such DC-specific surface structure is identified, an antibody thereto may be used alone in this stage of the protocol. At present, however, the use of antibodies in the second category of DC-distinguishing antibodies — antibodies that affirmatively bind to dendritic cell surface structures — obligates the additional use of DC-distinguishing antibodies from the first category, i.e., those that identify non-dendritic lineages.

Conversely, antibodies from the first category of DC-distinguishing antibodies — those that

preferentially bind non-dendritic cells - cannot at present be used without at least one antibody from the second category. Basophils are $\text{lin}^- \text{CD123}^{\text{high}} \text{CD11c}^+$ but HLA-DR $^-$; when antibodies that preferentially bind non-dendritic cells (category 1) are used in the assay, an anti-HLA-DR antibody must also be used.

If dendritic cell-distinguishing antibodies of both first and second category are used, the antibodies in the two categories are preferentially labeled with fluorophores that are flow cytometrically distinguishable.

The surface staining step may optionally also include, as a second broad class, antibodies that distinguish known dendritic cell subsets. Thus, antibodies specific for CD11c or CD123 prove particularly useful, as these antigens are known to define mutually exclusive peripheral blood DC subsets. The fluorophore used should be flow cytometrically distinguishable. Thus, where antibodies used later in the assay for intracellular staining are labeled with phycoerythrin (PE), a typical surface staining scheme would include, e.g., $\text{lin} 1$ FITC, HLA-DR PerCP, and CD11c APC (in this nomenclature, the antibody is identified by its specificity, followed by the fluorophore).

After surface staining, the red cells in the sample are lysed and the nucleated cells then permeabilized. These two steps may be accomplished using commercially available reagents, such as FACS[®] Permeabilizing Solution and FACS[®] Lysing solution (BDIS catalogue numbers 340457 and 349202, respectively), according to the manufacturer's instructions.

Following permeabilization, the cells are stained intracellularly using fluorophore-conjugated

antibodies that are specific for cytokines. The fluorophore conjugated to the cytokine-specific antibodies is preferentially distinguishable in a flow cytometric assay from any of those used for surface staining.

After intracellular staining, the cells are washed and then analyzed using a flow cytometer, preferably one capable of simultaneous excitation and detection of multiple fluorophores.

FIG. 1 does not schematize the assay for detecting changes in the surface expression of dendritic cell activation markers, which differs in some respects from that used to detect changes in cytokine expression.

In such an assay, activation of dendritic cells in whole blood is performed in the absence of secretion inhibitor, such as Brefeldin A. This precludes the concurrent measurement, in any such sample, of intracellular cytokine expression.

After incubation in the presence of activator, the surface of the cells is stained with fluorophore-conjugated antibodies. In this step, a plurality of dendritic-cell distinguishing antibodies is used, optionally with antibodies that distinguish known dendritic cell subsets, as above-described.

In addition, however, a third class of surface-staining antibodies is used. These are antibodies that recognize surface structures, typically proteins, the expression of which is altered by the prior incubation with dendritic cell activator. For example, activation of peripheral blood dendritic cells is known to cause upregulation of the T cell costimulatory molecules CD80 (B7.1), CD86 (B7.2) and HLA-DQ. Olweus et al., *Proc. Natl. Acad. Sci. USA*

94:12551-12556 (1997). Thus, the surface staining step, as desired, may include antibodies specific for one or more of these antigens. Recent reports identify CD83 and CMRF-44 as cell surface markers that are
5 expressed at high levels on activated or cultured DCs from blood and lymphoid tissue; antibodies specific for these markers may also advantageously be used. Antibodies of this class, if used, are typically conjugated to a fluorophore that is flow cytometrically
10 distinguishable from the antibodies described above. Thus, a typical surface staining scheme would include, e.g., lin 1 FITC, HLA-DR PerCP, CD11c APC, and an antibody specific to a DC surface activation antigen labeled with PE.

15 After surface staining, the red cells in the sample are lysed and the cells are washed and then analyzed using a flow cytometer, preferably one capable of simultaneous excitation and detection of multiple fluorophores.

20 As further elaborated in the experimental examples hereinbelow and FIGS. 2 - 8, whole blood samples from normal volunteers were assayed for dendritic cell function. Preparations were activated with either lipopolysaccharide ("LPS"), phorbol 12-myristate 13 acetate ("PMA") plus ionomycin ("I")
25 (together, "PMA + I") or CD40-crosslinking, each for 4 hours at 37°C. Substances attempted as activators that elicited no cytokine production - PHA, CD2/2R (BDIS Cat. No. 340366), SEB (staph enterotoxin B), CMV, and
30 crosslinking of CD49d - are not reported. CD40 crosslinking effected changes in surface antigen expression but failed to elicit cytokine production.

Table 1 lists the cytokines that were assayed in one or more of the experiments, further classified

according to the DC activator used in the experiment. A plus ("+") indicates that expression of the respective cytokine was assessed in one or more experiments; a minus ("-") indicates that expression of the respective cytokine was not assessed. Table 1 does not report the level of expression, which follows in Table 2.

10

15

20

	LPS	PMA + I	CD40 crosslinking
IL-1 α	+	+	+
IL-1 β	+	+	+
IL-1RA	+	+	+
IL-2	+	+	-
IL-4	+	+	-
IL-6	+	+	+
IL-8	+	+	+
IL-10	+	+	+
IL-12	+	+	+
IL-13	+	+	+
TNF α	+	+	+
IFN- γ	+	-	-

Table 2 presents the functional responses of the CD11c+ and CD11c- dendritic cell subsets, stimulated with either LPS or PMA+I, and assayed in whole blood. Cytokine expression is measured as mean fluorescence intensity (MFI); the change in surface molecule expression is measured as a ratio of mean

fluorescence intensities (MFI) of activated versus control sample.

Table 2 DC Functional Responses				
	LPS CD11c ⁺	PMA+I CD11c ⁺	LPS CD11c ⁻	PMA+I CD11c ⁻
5	IL-1 α	-	-/+	-
	IL-1 β	+	+	-
	IL-1RA	+	+	-
	IL-2	-	-	-
	IL-4	-	-	-
10	IL-6	+	-	-
	IL-8	+	+	-
	IL-10	-	-	-
	IL-12	-/+	-	-
	IL-13	-	-	-
15	TNF α	+	+	-
	IFN- γ	-	N.A.	N.A.
	CD25	+	-/+	+
	CD40	+	+	+
	CD80	+	+	-/+
20	CD86	+	+	-/+
	HLA-DQ	+	+	-
	HLA-DR	+	+	-/+

The results of these experiments, quite surprisingly, demonstrated the CD11c⁻CD123⁺ subset failed to produce any of the tested cytokines, no matter which DC activator was used. When assayed for

changes in surface antigen expression, this subset demonstrated clear upregulation of CD25 expression upon PMA+I activation; upregulation of CD25 was the only distinct response observed in CD11c⁻CD123⁺ DCs for all
5 investigated stimuli.

In striking contrast, the CD11c⁺CD123^{low} DCs showed easily measured changes in cytokine expression when stimulated with LPS or PMA+I.

With LPS stimulation, CD11c⁺ cells produced
10 high levels of TNF α and IL-1 β , lower levels of IL-6, IL-1RA and IL-8, and trace levels of IL-12 and IL-1a. The response to LPS is surprising: the CD11c⁺CD123^{low} DC are CD14⁻, and CD14 is the principal LPS receptor. It seems likely that LPS acts additionally through a
15 second receptor, perhaps CD11c itself. Consistent with that hypothesis, the CD11c⁻ DCs, which lack both CD14 and CD11c, fail to respond to LPS stimulation with increased intracellular cytokine expression. Further consistent with this hypothesis, CD14⁺CD11c⁺ monocytes
20 respond to LPS stimulation much more potently than do CD14⁻CD11c⁺ DCs (FIG. 7A), without showing significantly increased response to PMA+I (FIG. 7B).

FIG. 8 further shows the kinetics of the response of CD11c⁺ DCs to LPS. TNF α is produced first,
25 followed by IL-1 β and IL-6. CD80 is intensely upregulated after about 4 hours of activation.

With PMA+I activation, CD11c⁺ cells produced IL-8 and IL-1 β , lower but significant levels of IL-1RA and TNF α , trace amounts of IL-1 α , and no detectable
30 IL-6.

Thus, differences in the cytokine responses of the CD11c⁺ DC subset to various activators were readily observed. Principal among these differences is the expression of IL-6 uniquely when stimulated with

LPS, and the altered relative expression of IL-8 and TNF α .

The activation of CD11c⁺ DCs in whole blood also led to an increased expression of accessory
5 molecules. LPS activation triggered upregulation of CD25, CD40, CD80, CD86, HLA-DR and HLA-DQ. The T cell co-stimulatory molecules, in particular CD80, gave the strongest signal. PMA+I led to an upregulation of CD86, CD80, HLA-DQ and HLA-DR. Minimal increase of
10 CD25 and CD40 were observed.

Activation via crosslinking of CD40 resulted in increased levels of CD86, CD80, and minimal upregulation of HLA-DR.

These data, as further detailed in the
15 experimental examples that follow, demonstrate that peripheral blood DC subsets may readily be distinguished in whole blood by their differential production of cytokines and/or cell surface proteins in response to DC activators.

20 Furthermore, because the dendritic cells that were observed to respond to DC activators fall into a subset (CD11c⁺) known to be more potent in T cell activation than is the subset (CD11c⁻) showing no such response, the data further demonstrate that the
25 parameters measured in the method of the present invention – cytokine production and upregulation of surface activation antigens – directly correlate with DC function.

The ease with which the present invention
30 permits measurement of DC function in whole blood, without prior DC purification, was unexpected, because the low frequency of DCs in blood, coupled with the tendency of activated DCs to adhere to equipment, had

earlier suggested that too few events could be assayed in a blood sample of clinically-relevant size.

The ability to measure DC function in whole blood, without prior DC purification, offers
5 significant advantages.

From a procedural standpoint, the methods of the present invention eliminate the cell loss attendant upon all DC purification schemes, increasing sensitivity and reducing possible systematic bias.
10 Additionally, the minimal perturbation effected by the methods of the present invention reduces the chance for phenotypic changes resulting from experimental intervention. And as a flow cytometric assay, the methods of the present invention permit DC function to
15 be assessed on a cell-by-cell, rather than bulk, basis, permitting fine discrimination.

From the standpoint of the data made newly available by this invention, the methods of the present invention permit, for the first time, the ready and
20 rapid assessment of DC function in whole blood.

As applied to human patients, the methods of the present invention thus permit the measurement of DC function to be added to the existing roster of immune function assays, and will find utility in clinical
25 situations in which such existing immune function tests are presently used. For example, the methods of the present invention may advantageously be used, alone or in conjunction with flow cytometric quantitation of CD4⁺ T lymphocyte levels, in the clinical staging of AIDS
30 progression. The methods of the present invention may also be used, alone or in conjunction with existing assays, in the assessment of immune function in congenital, rather than acquired, immunodeficiency syndromes, and in the assessment of immune competence

following therapeutic immunosuppression or immunoablation. At the other end of the clinical spectrum, the methods of the present invention will also profitably find use, alone or in conjunction with
5 existing assays, in the clinical assessment of various forms of immune hypersensitivity, allergies, or in the clinical assessment of autoimmune diseases such as multiple sclerosis, rheumatoid arthritis, sarcoidosis, or the like.

10 By permitting the study of DC function in whole blood, the methods of the present invention also permit the ready evaluation of the effects that agents circulating in the blood may have on DC function. In particular, the assay permits the assessment of the
15 specific effects on DC function of pharmaceutical agents that either intentionally or fortuitously affect DC function.

Thus, as applied to the measurement of DC function in experimental mammals, whether outbred,
20 inbred, or transgenic, the methods of the present invention allow pharmaceutical agents to be tested for their *in vivo* effects on DC function, permitting the selection of agents that desirably demonstrate immunomodulatory effects, or the selection of agents
25 that specifically lack such effects.

As applied to human subjects, the methods of the present invention permit the ready assessment of the intentional or fortuitous effects on DC function of drugs that circulate in the patient's blood, as a
30 complement to existing immune function assays. For example, the methods of the present invention may be used to assist the monitoring and titration of immunosuppressive agents. The methods prove particularly useful in the monitoring and titration of

immunosuppressive agents that abrogate, downmodulate, or otherwise interfere with the function of cytokines, chemokines, or growth factors. Conversely, the methods of the present invention also prove particularly useful
5 in the monitoring of the effects of affirmative cytokine therapy, such as therapies involving the administration of interferons in the treatment of multiple sclerosis, the administration of growth factors after myeloablation, or the like.

10 The methods of the present invention may also be used to monitor immunomodulatory side effects of agents given to effect unrelated clinical goals.

 The methods of the present invention are particularly well suited to the experimental and
15 clinical assessment of therapies involving DC cells themselves. Thus, the methods of the present invention find use in the design, assessment, and monitoring of therapies in which autologous dendritic cells are administered after *in vitro* manipulation, therapies in
20 which dendritic cells are targeted for ablation, either *in vitro*, to facilitate transplantation, or *in vivo*, to effect immunosuppression or induction of tolerance, or therapies in which dendritic cells are targeted to increase global or specific immune function.

25 It will be understood that the dendritic cells that are found circulating in the peripheral blood at the moment that blood is drawn — those assayed in the methods of the present invention — are drawn from a temporal window in the maturation of one or more
30 cell lineages. That is, for each discrete lineage, the cells circulate preferentially during particular phases in the maturational process. Nothing, however, in the methods of the present invention is limited to particular phases in DC maturation. Thus, the methods

may equally be applied to CD34⁺ committed DC precursors that circulate spontaneously, or to CD34⁺ DC precursors that are mobilized by pharmacological intervention or the like.

5 The invention may be better understood by reference to the following examples, which are offered by way of illustration and not by way of limitation.

EXAMPLE 1 **Materials**

10 Unless otherwise specified, the following reagents were used in the experiments presented herein. For convenience, antibodies are identified by their specificities and conjugated fluorophore. Fluorophores are phycoerythrin (PE), peridinin chlorophyll protein
15 (PerCP), allophycocyanin (APC), fluorescein isothiocyanate (FITC). Thus, an antibody labeled with phycoerythrin (PE) that is specific for TNF α is denominated "TNF α PE".

Antibodies

20 The following antibodies were obtained from Becton-Dickinson Immunocytometry Systems (BDIS, San Jose, CA): TNF α PE; IL-1 α PE; IL-1RA PE; IL-1 β PE; IL-2 PE; IL-4 PE; IL-6 PE; IL-8 PE; IL-13 PE; IFN- γ PE; CD11c APC (5 μ L/test) (0.125 μ g/test); HLA-DR PerCP (10
25 μ L/test) (0.125 μ g/test); lin 1 FITC (research lot KW98/07 1.1) (20 μ L/test); CD40 PE (unconjugated mAb obtained from DNAX Research Institute, Palo Alto, CA; custom conjugated to PE at BDIS, BDIS research conjugate PC#931) (10 μ L/test) (0.125 μ g/test); CD80 PE

(20 μ L/test); CD25 PE (20 μ L/test); HLA-DQ PE (unconjugated mAb obtained from BDIS; custom-conjugated to PE at BDIS, BDIS research conjugate PC#1284) (0.5 μ g/test); IgG2a PE (cat. no. 340459, 25 μ g/mL); IgG1 PE (cat. no. 340013, 50 μ g/mL). CD83 PE (BDIS research conjugate, PC#1520; 50 μ g/mL, used at 20 μ L per 300 μ L whole blood) is also available commercially from Pharmingen (catalogue number 36935X).

The lin 1 FITC lineage cocktail is also available commercially (BDIS, catalogue number 340546), and contains a titrated mixture of antibodies specific for CD3, CD14, CD16, CD19, CD20, and CD56, all labeled with FITC. In combination, the antibodies stain lymphocytes, monocytes, eosinophils, and neutrophils.

The following antibodies were obtained from PharMingen (San Diego, CA): IL-10 PE (IgG2a) (0.1 μ g/test); IL-12 PE (IgG1) (0.1 μ g/test); CD86 PE (clone IT2.2; Cat #33435X, IgG2b) (10 μ L/test).

20 Lysing and permeabilizing agents

FACS® Permeabilizing Solution and FACS® Lysing Solution were obtained as 10X stock solutions from BDIS (catalogue numbers 340457 and 349202, respectively), and were diluted and used in accordance with the package insert.

Dendritic cell activators

Chemical activators were obtained from Sigma Chemical Company, St. Louis, MO. Lipopolysaccharide ("LPS") (Sigma catalogue number L2654) was made 0.5 mg/mL in DMSO and stored at -20°C. Ionomycin ("I") (Catalogue number I-0634), was made 0.5 mg/mL in

ethanol and stored at -20°C. Phorbol 12-myristate 13 acetate ("PMA") (Catalogue number P-8139) was made 0.1 mg/mL in DMSO and stored at -20°C.

CD40 crosslinking was performed using
5 polystyrene beads (0.84 µm, Baxter) coated with CD40 antibody (PharMingen, San Diego).

Secretion inhibitor

Brefeldin A ("BFA") (catalogue number B-7651) was made 5 mg/mL in DMSO, and stored at -20°C.

10 Wash buffer

Wash buffer consisted of phosphate-buffered saline ("DPBS") (obtained as a 10X stock solution from GibCoBRL (Grand Island, NY), then diluted with deionized water to 1X), containing 0.5% fetal calf
15 serum (Sigma, St. Louis, MO) (fetal calf serum added after dilution of 10X PBS stock to 1X).

EXAMPLE 2

20 Protocols for Whole Blood Flow Cytometric Dendritic Cell Immune Function Assay

Unless otherwise specified, the following protocols were used in the experiments presented herein.

Dendritic cell activation

25 Venous blood of normal donors was collected in sodium heparin VACUTAINER® tubes. For activation with LPS, the blood was stimulated with 1 µg/mL LPS.

For activation with PMA+I, whole blood was first diluted 1:1 with RPMI medium (Biowhittaker, Watersville, MD). PMA was then added at 5 ng/mL and ionomycin at 1 µg/mL. For activation by CD40 crosslinking, 50 µL CD40-coated polystyrene beads was added to 1 mL whole blood. All samples were incubated for four hours at 37°C in a humidified incubator with 5% CO₂.

For detection of intracellular cytokines, activation, as above, was performed in the presence additionally of Brefeldin A (BFA) at 10 µg/mL. Control (resting) aliquots were incubated with BFA alone.

For detection of changes in surface antigen expression, samples were incubated with DC activator, as above, without the further addition of BFA. Control (resting) aliquots were incubated with neither BFA nor activator.

Immunofluorescence staining of intracellular cytokines

Prior to staining, PMA+I treated blood samples were reduced to half volume by centrifugation and removal of supernatant.

Cell preparation was done at room temperature (RT), and all incubation steps were performed in the dark.

For staining, 1 mL of sample (activated or resting blood control) was added to a cocktail of dendritic cell-distinguishing antibodies (20 µL lineage cocktail 1-FITC, 10 µL HLA-DR PerCP, 5 µL CD11c APC; reagent volumes per 50 µL blood) in a 50 mL polypropylene centrifuge tube. The blood was incubated in the presence of the fluorophore-conjugated antibodies for 15 min. After incubation, 40 mL FACS®

Lysing Solution was added and the tube incubated for a further 10 min. The cells were then collected by centrifugation for 10 min at 500 x g, and the pellet gently broken off for further processing. Next, 10 mL
5 FACS® Permeabilizing Solution was added and the cells were incubated for 10 min. The permeabilization reaction was stopped by addition of 40 mL of buffer (DPBS 1X, 0.5% fetal calf serum). The permeabilized cells were pelleted for 10 min at 500 x g and
10 resuspended in the supernatant remaining in the tube after decanting (approximate volume 500 µL).

An aliquot of 50 µL of the extracellularly-stained, lysed and permeabilized cells (sufficient for one test) was added to a polypropylene staining tube
15 and incubated for 30 min in the presence of the cytokine-specific mAb (see Materials, above). The samples were then washed with buffer, resuspended in 250 µL buffer, and subjected to flow cytometric data acquisition as soon thereafter as possible. If flow
20 cytometric data acquisition was delayed, the samples were kept at 4°C for up to one hour.

Depending upon yield, a 1 mL sample of whole blood yielded about 7 to 12 tests for cytokine expression determination.

25 Immunofluorescent staining of surface antigens

Prior to staining, PMA+I treated blood samples were reduced to half volume by centrifugation and removal of supernatant.

Cell preparation was done at room temperature
30 (RT) and all incubation steps were performed in the dark.

For staining, 150 μ L of sample (activated or resting blood control) was added to a cocktail of monoclonal antibodies in a staining tube. The cocktail included a plurality of dendritic-cell distinguishing antibodies (20 μ L lin 1 FITC cocktail, 10 μ L HLA-DR PerCP, 5 μ L CD11cAPC; reagent volumes per 100 μ L blood) and one of the following PE-conjugated antibodies specific for DC surface activation antigens (20 μ L CD25 PE, 0.125 μ g CD40 PE, 20 μ L CD80 PE, 10 μ L CD86 PE, 10 0.5 μ g HLA-DQ PE). Blood and mAbs were incubated for 15 min at RT in the dark.

After incubation, 3 mL of FACS® Lysing Solution was added and the tube incubated for 10 min at RT. The lysed cells were centrifuged for 5 min at 15 500 x g and subsequently washed with 3 mL buffer (DPBS 1X, 0.5% fetal calf serum). The cell pellet was resuspended in 250 μ L buffer and immediately acquired on a flow cytometer. If data acquisition was delayed, the cells were maintained at 4°C for up to one hour.

20

Flow Cytometric Analysis

The samples as described above were acquired on a FACSCalibur™ dual laser flow cytometer (BDIS, San Jose, CA). The instrument was set up using automated 25 FACSComp™ 4.0 software and 4-color Calibrite™ beads (BDIS, San Jose, CA). Events were acquired on a FSC threshold. To reduce the size of the listmode data files, the acquisition used a live gate on HLA-DR positive events in a lin 1 FITC/HLA-DR PerCP two- 30 parameter distribution.

EXAMPLE 3

Detection of CD11c⁺ DC Cytokine Response in Whole Blood

Whole blood samples were drawn from healthy volunteers and activated with either LPS or PMA+I, both
5 in the presence of Brefeldin A, according to the procedures described in Examples 1 and 2. The results are shown, respectively, in FIGS. 2 and 3.

FIGS. 2A - 2C show the surface immunophenotypic characteristics of peripheral blood DC
10 from a single LPS-activated whole blood sample. CD11c⁺ dendritic cells are painted green, CD11c⁻ DC are painted red, and nondendritic cells appear gray. The colors are arbitrarily chosen for purposes of display, and bear no relationship to the fluorophores used for
15 analysis. FIG. 2A demonstrates that both dendritic cell subsets are lin¹ FITC^{dim} and HLA-DR^{bright}, in agreement with O'Doherty et al., *Immunology* 82:487-493 (1994); Olweus et al., *Proc. Natl. Acad. Sci. USA* 94(23): 12551-12556 (1997), with FIG. 2B further
20 demonstrating that the two subsets have similar side scatter and forward scatter properties. FIG. 2C shows discrimination of the two subsets based on differential levels of CD11c expression.

FIGS. 2D - 2J show the result of assays for
25 expression of IL-1RA (FIG. 2D), TNF α (FIG. 2E), IL-6 (FIG. 2F), IL-8 (FIG. 2G), IL-12 (FIG. 2H), IL-1 α (FIG. 2I). FIG. 2J shows results using an isotype-matched PE-conjugated negative control antibody.

FIGS. 2D - 2J demonstrate that the CD11c⁻
30 (CD123⁺) subset (red) is unresponsive to LPS stimulation, at least as evidenced by the absence of detectable cytokine production. Although not shown directly on these figures, the cytokine levels measured

in the LPS-activated CD11c⁻ DC are indistinguishable from those produced in the absence of activator; as shown in FIG. 4, neither CD11c⁻ nor CD11c⁺ subset produces detectable levels of cytokine in the absence of DC activators.

In contrast, the CD11c⁺ population shows much higher levels of cytokine production, with high levels of TNF α and IL-1 β , lower levels of IL-6, IL-1RA and IL-8, and trace levels of IL-12 and IL-1a.

FIGS. 3A - 3C show the surface immunophenotypic characteristics of peripheral blood DC from a single whole blood sample activated with PMA+I. CD11c⁺ dendritic cells are painted green, CD11c⁻ DC are painted red, and nondendritic cells appear gray. The colors are arbitrarily chosen for purposes of display, and bear no relationship to the fluorophores used for analysis.

FIG. 3A demonstrates that both dendritic cell subsets are lin 1 FITC^{dim} and HLA-DR^{bright}, with FIG. 3B further demonstrating that the two subsets have similar side scatter and forward scatter properties. FIG. 3C shows discrimination of the two subsets based on differential levels of CD11c expression.

FIGS. 3D - 3I show the result of assays for expression of TNF α (FIG. 3D), IL-1 α (FIG. 3E), IL-1 β (FIG. 3F), IL-1RA (FIG. 3G), and IL-8 (FIG. 3H). FIG. 3I shows results using an isotype-matched PE-conjugated negative control antibody.

FIGS. 3D - 3I demonstrate that the CD11c⁻ (CD123⁺) subset (red) is unresponsive to PMA + I stimulation, at least as evidenced by the absence of detectable cytokine production. Although not shown directly on these figures, the cytokine levels measured in the LPS-activated CD11c⁻ DC are indistinguishable

from those produced in the absence of activator
(compare to FIG. 4).

In contrast, the CD11c⁺ population shows much
higher levels of cytokine production, with demonstrable
production of IL-1 β , IL-1RA, TNF α , and IL-8. FIG. 3E
demonstrates that CD11c⁺ cells produced trace amounts of
IL-1 α .

The cytokine expression of monocytes was
evaluated in some of the same samples. Monocytes were
identified based on their scatter characteristics,
their bright lin 1 FITC, anti-HLA-DR PerCP and CD11c
APC staining. As can be seen, in LPS+BFA stimulated
samples, monocytes express cytokines at higher levels
than do CD11c⁺ DCs. Upon PMA+I+BFA activation, the
cytokine secretion of CD11c⁺ DCs and monocytes is
equivalent, but much less compared to LPS activated
samples. Intracellular protein secretion is evaluated
as PE mean fluorescence intensity (MFI).

EXAMPLE 4 Detection of CD11c⁺ DC Surface Antigen Expression in Whole Blood

Whole blood samples were drawn from healthy
volunteers and activated with LPS, PMA+I, or CD40
crosslinking, in the absence of Brefeldin A, according
to the procedures described in Examples 1 and 2. The
results are shown in FIG. 6.

As demonstrated in the histograms, the CD11c⁺
subset demonstrated clear upregulation of CD25
expression upon PMA+I activation; upregulation of CD25
was the only distinct response observed in CD11c⁺
subset. In contrast, the CD11c⁺ subset showed
upregulation of CD25, CD40, CD80, CD86, HLA-DR and HLA-
DQ upon LPS activation. The T cell co-stimulatory

molecules, in particular CD80, gave the strongest signal. PMA+I led to an upregulation in CD11c⁺ cells of CD86, CD80, HLA-DQ and HLA-DR. Minimal increase of CD25 and CD40 were observed. Activation via
5 crosslinking of CD40 resulted in increased levels of CD86, CD80, and minimal upregulation of HLA-DR.

FIG. 5 highlights the differences in the TNF α , IL-8, CD80 and CD86 responses of the CD11c⁺ peripheral blood DC subset during activation with PMA+I
10 versus LPS. Two donors are displayed. The expression pattern of cytokines and co-stimulatory molecules varies between different stimuli and is consistent between donors. IL-8 and CD86 are the dominant signal in PMA+I stimulation. In LPS activated samples, TNF α
15 and CD80 are produced to a greater extent.

EXAMPLE 5 Cytokine Kinetic Assay

Whole blood was stimulated with 1 μ g/ml LPS and incubated in polypropylene tubes (12 x 75 mm) at
20 37°C in a humidified atmosphere containing 5% CO₂. Aliquots were processed, essentially as set forth in Examples 1 and 2 above, every hour from 0 to 8 hrs post stimulation. One hour prior to each sample harvest, BFA was added at 10 μ g/mL. Results are shown in FIG.
25 8.

EXAMPLE 6

Modified Protocols for Whole Blood Flow Cytometric
Dendritic Cell Immune Function Assay

This Example presents protocols modified
5 slightly from those set forth in Examples 1 and 2; the
modified protocols are faster, permit data to be
acquired from a greater number of cells, thus improving
statistics, and contain fewer aliquoting steps.

DC Surface Antigen Expression in Whole Blood

10 Response of CD11c⁺ DC to LPS stimulation,
reported as the mean fluorescence intensity (mean MFI)
of surface expression of the co-stimulatory molecule
CD80 or the DC activation marker CD83, is measured as
follows.

15 Materials

- lipopolysaccharide (LPS), by Sigma,
Cat#L2654, 0.5 mg/mL in DMSO, stored at
-70°C, working solution is a [1:100] dilution
in PBS (1x)
- 20 • DMSO [1:100] in PBS (1x) (stimulus solvent)
- CD11c APC, by BDIS, Cat#340544, use 15 µL/300
µL whole blood (WB)
- lineage cocktail 1 (lin 1) FITC, by BDIS,
Cat#340546, 60 µL/300 µL WB
- 25 • Anti-human HLA-DR PerCP, by BDIS, Cat#347364,
use 30 µL/300 µL WB

- CD80 PE, by BDIS, Cat#340512, use 60 μ L/300 μ L WB
- CD83 PE, by BDIS (custom conjugate), PC#1520, 50 μ g/mL, use 20 μ L/300 μ L WB
- 5 • FACS™ lysing Solution 10x, by BDIS, Cat #349202, working solution is a [1:10] dilution in deionized water
- wash buffer PBS (1x), 0.5% BSA
- PBS (1x)

10 Stimulation protocol

- Add 20 μ L of LPS working solution to 1 mL of whole human blood (activated control).
- Add 20 μ L of DMSO:PBS1x [1:100] to 1 mL of whole human blood (resting control).
- 15 • Incubate tubes for 4 hrs @ 37°C and 5% CO₂.

Staining protocol

Perform all steps in 12x75 polypropylene tubes.

- add 60 μ L lin 1 FITC, 30 μ L anti-Hu HLA-DR PerCP, 15 μ L CD11c APC and 20 μ L CD83 PE or
20 60 μ L CD80 PE
- add 300 μ L whole blood (WB), obtained from the stimulation protocol, vortex

- incubate for 15 min @ RT in the dark
- add 3 mL FACSlysing Solution (1x), cap tubes and vortex
- incubate for 10 min @ RT in the dark
- 5 • pellet (7 min @ 1500 rpm, Sorvall benchtop centrifuge)
- aspirate supernatant and vortex to break off pellet
- add 2 mL wash buffer and vortex
- 10 • pellet (7 min @ 1500 rpm, Sorvall benchtop centrifuge)
- aspirate supernatant
- resuspend stained cells in ~200 μ L wash buffer
- 15 • acquire samples (i.e., perform flow cytometric analysis) in \leq 1 hr, store samples prior to acquisition @ 4°C

20 Note: It is necessary to transfer samples into 12x75 polystyrene tubes for the acquisition; polypropylene tubes do not commonly fit the sample injection system of the flow cytometer. Perform the transfer right before acquisition.

Instrument & Setup & Acquisition Criteria

- FACSCalibur Flow Cytometer
 - 4-Color Lyse/ Wash (LW) FACSCComp software setting
 - Use a FSC threshold
- 5 • Acquire the FCS file using anti-human HLA-DR PerCP positive / lin 1 FITC dim events. Use the entire cell suspension for acquisition (about 4,000 -7,000 events).

Cell Identification

- 10 Scatter low / lin 1 FITC low / anti-human HLA-DR PerCP high / CD11c APC high

DC Cytokine Expression in Whole Blood

- Response of CD11c⁺ DC to LPS stimulation, reported as mean fluorescence intensity (mean MFI) of
- 15 TNF α expression in the presence of the Golgi-transport disrupting agent brefeldin A (BFA), is measured as follows. The cells are activated for 2 hrs.

Materials

- 20 • lipopolysaccharide (LPS), by Sigma, Cat#L2654, 0.5 mg/mL in DMSO, stored at -70°C, working solution is a [1:100] dilution in PBS (1x)
- DMSO [1:100] in PBS (1x) (stimulus solvent)

- Brefeldin A (BFA), by Sigma, Cat# B-7651, stock solution 50 mg/mL in DMSO, stored at -20°C, working solution [1:100] in PBS(1x), use 20 µL/1 mL WB
- 5 • CD11c APC, by BDIS, Cat#340544, use 30 µL/300 µL whole blood (WB)
- lineage cocktail 1 (lin 1) FITC, by BDIS, Cat#340546, 120 µL/300 µL WB
- 10 • Anti-human HLA-DR PerCP, by BDIS, Cat#347364, use 60 µL/300 µL WB
- Anti-human TNFα PE, by BDIS, Cat#340512, use 20 µL/test
- 15 • FACS™ Permeabilizing Solution 10x, by BDIS, Cat#340457, working solution is [1:10] dilution in deionized water
- FACS™ lysing Solution 10x, by BDIS, Cat#349202, working solution is a [1:10] dilution in deionized water
- wash buffer PBS(1x), 0.5% BSA
- 20 • PBS(1x)

Stimulation protocol

- Add 20 µL of LPS working solution to 1 mL of whole human blood (activated control).

- Add 20 μ L of DMSO:PBS1x [1:100] to 1 mL of whole human blood (resting control).
 - Then add 20 μ L of BFA working solution to both control tubes.
- 5 • Incubate tubes for 2 hrs @ 37°C and 5% CO₂.

Staining protocol (intracellular)

- Perform all steps in 12x75 polypropylene tubes.
- 10 • add 120 μ L lin 1 FITC, 60 μ L anti-Human HLA-DR PerCP, 30 μ L CD11c APC
- add 300 μ L whole blood (WB), obtained from the stimulation protocol, vortex
- incubate for 15 min @ RT in the dark
- 15 • add 3 mL FACS lysing Solution 1x, cap tubes and vortex
- incubate for 10 min @ RT in the dark
- pellet (7 min @ 1500 rpm, Sorvall benchtop centrifuge)
- 20 • aspirate supernatant and vortex gently to break off pellet
- resuspend in 1 mL FACS permeabilizing Solution 1x

- incubate @ RT for 10 min in the dark
- add 2.5 mL wash buffer, cap tubes and vortex gently
- 5 • pellet (10 min @ 1500 rpm, Sorvall benchtop centrifuge)
- decant supernatant
- resuspend permeabilized cell pellet gently in remaining supernatant
- 10 • add 20 μ L anti-human TNF α PE reagent for intracellular staining
- incubate for 30 min @ RT in the dark
- add 2 mL wash buffer and gently vortex
- pellet (10 min @ 1500 rpm, Sorvall benchtop centrifuge)
- 15 • decant supernatant
- resuspend stained and permeabilized cells in ~200 μ L wash buffer
- acquire sample in \leq 1 hr, store samples prior to acquisition @ 4°C
- 20 Note: It is necessary to transfer samples into 12x75 polystyrene tubes for the acquisition; polypropylene tubes do not commonly fit the sample

injection system of the flow cytometer. Perform the transfer right before acquisition.

Instrument & Setup & Acquisition Criteria

- FACSCalibur Flow Cytometer.
- 5 • 4-Color Lyse/No Wash (LNW) setting in FACSComp software.
- Use a FSC threshold.
- 10 • Acquire the FCS file using anti-human HLA-DR PerCP positive / lin 1 FITC dim events. Use the entire cell suspension for acquisition (4,000 -7,000 events).

Cell Identification

- 15 Scatter low / lin 1 FITC low / Anti-human HLA-DR PerCP high / CD11c APC high.

20 All patents, patent publications, and other published references mentioned herein are hereby incorporated by reference in their entirety as if each had been individually and specifically incorporated by reference herein.

25 While preferred illustrative embodiments of the present invention are described, it will be apparent to one skilled in the art that various changes and modifications may be made therein without departing from the invention, and it is intended in the appended

claims to cover all such changes and modifications which fall within the true spirit and scope of the invention.